

Drinkwell 360 replacement filters



DRINKWELL 360

Replacement Filters 3-Pack

Only for use with the Drinkwell® 360 Fountain

Filters Keep Water Clean and Fresh

- Bio Density Polyester Pad captures debris
- Granulated Carbon for extended action

PAC08-13069

40744Q

Only for use with the Drinkwell® 360 Fountain

www.petfountain.com

DRINKWELL 360

Replacement Filters 3-Pack

Filters Keep Water Clean and Fresh

- Bio Density Polyester Pad captures debris
- Granulated Carbon for extended action

THE ORIGINAL PET FOUNTAIN™

DRINKWELL 360

Replacement Filters 3-Pack

The constant circulation and filtration in the Drinkwell® Pet Fountain removes bad tastes and odors, making the water fresher and more appealing to your pet.

Veterinary Ventures, Inc.
PO BOX 19717
Berno, NY 89511

www.petfountain.com

DRINKWELL 360

Replacement Filters 3-Pack

Only for use with the Drinkwell® 360 Fountain

Instructions:

Rinse the filter thoroughly under cool water before placing into the filter housing. This will remove any loose charcoal dust, which is harmless. A small amount of charcoal dust may shed for the life of the filter. This is normal and not harmful to your pet.

Remove the Flow Control Cap and sprouting. Then remove the Filter Housing Cap and insert the Bio-Density® 360 Filter. Push in a pull tab which should be placed in the filter housing with the tab facing upwards. In every placement area, removal. Next, screw on the filter housing lid to secure the filter in place. See product manual for detailed instructions or visit us at: www.petfountain.com

Filters should be changed every 2 weeks, depending on the size and number of pets using the fountain.