



Replacing the tubing

After months of using your fountain, you may want to change the tubing out. In the long run, the most economical way to do this is to buy a 100 foot role (Hardware World offers it for \$19.72, \$25.73 w/ shipping). Below is a Google search; you are looking for Samar beverage grade tubing I.D. (inner diameter) ¼ inch, O.D. (outer diameter) 3/8 inch. It's important to go with Samar (see below) because they are a US manufacturer of FDA-compliant beverage grade tubing (for the safety of our kitties). The tubing at your local hardware store may be fine, too, but we prefer to go with certified safety.

The vinyl PVC tubing we use:



Samar's versatile Clear Vinyl Tubing can be used for beverage dispensing, potable water, laboratory tubing, and numerous low pressure applications. It is light weight, yet tough and abrasion resistant. The CVT's PVC compound is in compliance with FDA CFR TITLE 21 requirements, and meets USDA, 3A, and NSF51 criteria.

All fountains use Samar 7003P from pump to filter ([Samar type 7003P 3/8" OD](#))

The newer faucet spout fountains use Samar 7004P (7/16" OD) ([Samar type 7004P 7/16" OD](#)) or 7005P (1/2" OD) ([Samar type 7005P 1/2" OD](#)) from filter to spout.

All other fountains use ([Samar type 7003P 3/8" OD](#)) from filter to spout.

You can go to our website (<http://ebifountains.com/petpumps.html>) for a direct link to an order site:

Samar type 7003P: <http://www.sears.com/samar-company-inc-clear-vinyl-tubing-7003p-100ft-3-8-od-vinyl/p-SPM6582454701>

Samar type 7004P: <http://www.sears.com/samar-company-7-16x100-nontoxic-vinyl-tubing/p-SPM6583226301?prdNo=1&blockNo=1&blockType=G1>

Samar type 7005P: <http://www.sears.com/samar-company-1-2x100ft-nontoxic-vinyl-tubing-by-samar/p-SPM6582452301?prdNo=9&blockNo=9&blockType=G9>

or use Google to search for a better price:

Google search for:

3/8 od clear vinyl beverage tubing Samar 7003p

7/16 od clear vinyl beverage tubing Samar 7004p

Samar 3/8 x 1/2 pvc Tubing 7005p Vinyl